



P 2.3.4

Conversion of electrical energy

- P 2.3.4.1 Converting electrical into heat energy – measuring with the voltmeter and ammeter
- P 2.3.4.2 Converting electrical into heat energy – measuring with the joule and wattmeter
- P 2.3.4.3 Converting electrical into heat energy – measuring with CASSY

Converting electrical into heat energy - measuring with the joule and wattmeter (P 2.3.4.2)

Cat. No.	Description	P 2.3.4.3					
		P 2.3.4.1 (a)	P 2.3.4.1 (b)	P 2.3.4.1 (c)	P 2.3.4.2 (a)	P 2.3.4.2 (b)	P 2.3.4.2 (c)
384 20	Electric calorimeter attachment						
386 48	Dewar vessel calorimeter with base						
388 02	Copper-block calorimeter with heating coil				1	1	1
388 03	Aluminum-block calorimeter with heating				1	1	1
388 04	Large aluminum-block calorimeter with heating coil				1	1	1
388 06	Pair of connecting cables				1	1	1
521 35	Voltage source, 0 ... 12 V, e.g. Variable extra-low voltage transformer S				1	1	1
382 34	Thermometer, -10° to + 110 °C						
388 05	Thermometer for calorimeters				1		
666 190	Digital thermometer with 1 input					1	
666 193	Temperature sensor NiCr-Ni					1	
531 120	Voltmeter, AC, U < 12 V, e.g. Multimeter LDanalog 20						
531 130	Ammeter, I < 6 A, e.g. Multimeter LDanalog 30						
313 07	Stopclock I, 30s/15min						
664 103	Beaker, 250 ml, ss, hard glass						
665 755	Graduated cylinder, 250 ml: 2						
501 28	Connecting lead, Ø 2.5 mm², 50 cm, black				3	3	3
501 45	Pair of cables, 50 cm, red and blue				1	1	1
531 831	Joule and wattmeter				1	1	1
524 009	Mobile-CASSY				1	1	1
524 0673	NiCr-Ni Adapter S				1	1	1
529 676	NiCr-Ni temperature sensor 1.5 mm				1	1	1
524 010USB	Sensor-CASSY				1	1	1
524 200	CASSY- Lab				1	1	1
additionally required:							
1 PC with Windows 95/NT or higher							

Just like mechanical energy, electrical energy can also be converted into heat. We can use e.g. a calorimeter vessel with a wire winding to which a voltage is connected to demonstrate this fact. When a current flows through the wire, Joule heat is generated and heats the calorimeter.

The supplied electrical energy

$$W(t) = U \cdot I \cdot t$$

is determined in the first experiment by measuring the voltage U , the current I and the time t , and in the second experiment measured directly using the joule and wattmeter. This results in a change in the temperature of the calorimeter which corresponds to the specific heat capacity

$$Q(t) = m \cdot c \cdot (\vartheta(t) - \vartheta(0)),$$

c : specific heat capacity, m : mass, $\vartheta(t)$: temperature at time t

To confirm the equivalence

$$Q(t) = W(t)$$

the two quantities are plotted together in a diagram.